


NovAseptic® Valve

Manual Stainless Steel Actuator, Spring Closed

Catalogue Number NA##/170

Manual Stainless Steel Actuator, Spring Closed Configuration



Nominal Dimensions in mm (in.)

NA##	Size	A	B	C	D*	E	F	Compatible Valve Body	
								NA##	NU###
NA12	½"	61.0 (2.402)	50.0 (1.949)	47.0 (1.850)	TR40x2	M5	5.0 (0.197)	NA12	NU050

* Note: Non-standard thread, Actuator specific.

Specifications

Net Weight (approximate)

Actuator NA12

Weight kg (lb) 0.4 (0.88)

Material

	Cover Stainless Steel in Compliance with		Indicator Stainless Steel in Compliance with	Handle Stainless Steel in Compliance with
Material Code	EN 1.4305	ASTM® 303	SS 2346	SS 2346

Surface Roughness: External surface $R_a \leq 1.6\mu\text{m}$ (63µin)

Environment Temperature: 0 to 70 °C (32 to 158 °F)

Design Pressure: -1 to 10 bar(g) (-14.5 to 145 psi(g))

Note! The applied valve body and diaphragm may have different design temperature and/or pressure limits. The weakest component in the assembled product determines the maximum design temperature and pressure limits.

Labeling: Each Actuator is labeled for full LOT traceability according to Millipore QA routines.

Packaging: The Actuator is packaged in a closed box.

Quality Control: Millipore Quality Assurance System guarantees the control and traceability of the product.

Options: For non-standard NovAseptic Actuator Options, please contact Millipore for further information.

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office or visit us at www.millipore.com.

In the U.S. and Canada, call toll-free 1-800-MILLIPORE (1-800-645-5476). In the U.S., Canada and Puerto Rico, fax orders to 1-800-MILLIFX (1-800-645-5439).

MILLIPORE

Millipore and NovAseptic are registered trademarks of Millipore Corporation or an affiliated company.

ASTM is a registered trademark of the American Society for Testing and Materials.

Specifications subject to change without notice.

Lit. No. SP1119EN00 Rev. -, 09/2006. Replaces DS-063

© 2006 Millipore Corporation. All rights reserved.