

Specification Sheet


NovAseptic® Valve, Diaphragm, Fluoroelastomer

Catalogue Number:

NA#/24

Product Description

The NovAseptic® Valve is specifically designed for aseptic applications and complies with the most stringent cleanability and sterilization requirements. The focus on aseptic design is a significant feature for all valves. The NovAseptic® Valve is designed for minimum dead leg, complete drain ability, and high chemical resistance.



Catalogue No. Structure

NA# / 24
Diaphragm Model 24 = Fluoroelastomer Diaphragm

Nominal Dimensions in mm (in.)

	A	B	ØC	Compatible Valve NA#	Compatible Valve NU#
NA12/24	34.0 (1.339)	16.0 (0.63)	18.0 (0.709)	NA12	NU050
NA18/24	46.5 (1.831)	23.0 (0.906)	25.0 (0.984)	NA18	NU075
NA25/24	62.0 (2.441)	29.3 (1.154)	34.0 (1.339)	NA25	NU100
NA38/24	75.0 (2.953)	42.4 (1.669)	48.3 (1.902)	NA38	NU150
NA51/24	97.0 (3.819)	52.4 (2.063)	64.0 (2.520)	NA51	NU200

*Note: Non-standard thread, special to fit NovAseptic® Actuators

Merck Millipore is a division of 

Specifications

Net Weight (approximate)					
Diaphragm	NA12	NA18	NA25	NA38	NA51
Weight kg (lb)	0.01 (0.022)	0.02 (0.044)	0.05 (0.11)	0.10 (0.220)	0.16 (0.353)

Design Temperature	
	Short Term Use*
Max (Dry Heat)	200 °C (392 °F)
Max (Steam)	150 °C (302 °F)
Min	(-)20 °C ((-)4 °F)

Additional Information	
Material	The Diaphragm is manufactured from 100% Fluorocarbon rubber.
Surface Roughness	Smooth (hydrophobic)
Design Pressure	(-) 1 to 6 bar (g) (- 14.5 to 87.0 psi (g))
Labeling	Each Diaphragm is labeled for full LOT traceability according to QA routines.
Packaging	The Diaphragm is packaged in a closed box.
Quality Control	Our Quality Assurance System guarantees the control and traceability of the product.
In Compliance with	All Diaphragms are manufactured with materials in compliance with FDA regulations §177.2600 and supplied with statement. The Diaphragm meets USP Biological Tests for Plastics, Class VI

*Note: The applied valve body and actuator may have different design temperature and/or pressure limits. The weakest component in the assembled product determines the maximum design temperature and pressure limits.



www.merckmillipore.com

Merck Millipore and the M mark are trademarks of Merck KGaA, Darmstadt, Germany.
NovAseptic is a registered trademark of Merck KGaA, Darmstadt, Germany.
Specifications subject to change without notice.
Lit. No. SP1021EN00 Rev. C 03/2012 DP SBU-11-05530 Printed in the USA.
©2012 Merck Millipore Corporation, Billerica, MA 01821 USA. All rights reserved.

To Place an Order or Receive Technical Assistance

In Europe, please call Customer Service:

France: 0825 045 645

Germany: 01805 045 645

Italy: 848 845 645

Spain: 901 516 645 Option 1

Switzerland: 0848 645 645

United Kingdom: 0870 900 4645

For other countries across Europe,
please call: +44 (0) 115 943 0840

Or visit www.merckmillipore.com/offices

For Technical Service visit:

www.merckmillipore.com/techservice